

Rapid Identification And Confirmation Of Extended Spectrum Beta-Lactamase Utilizing *Brilliance*[™] ESBL Agar And Sensititre[™] ESBL Antimicrobial Susceptibility Panels

P. Stagg¹, D. Paisey¹, Bains, J¹, A. Appleton², C. Bastulli³. ¹Thermo Fisher Scientific, East Grinstead, UK, ²Thermo Fisher Scientific, Basingstoke, UK, ³Thermo Fisher Scientific, Cleveland Ohio, US

Overview

Purpose: This study was undertaken to determine if Sensititre[™] ESBL antimicrobial susceptibility panels could be used directly with *Brilliance*[™] ESBL agar for the rapid identification, confirmation and MIC determination of suspect organisms.

Methods: A collection of 66 clinical isolates were inoculated onto *Brilliance* ESBL Agar and Columbia Blood Agar (CBA). Post-incubation, Sensititre ESBL antimicrobial susceptibility testing panel inocula were prepared from colonies grown on both *Brilliance* ESBL Agar and CBA plates.


Results : The Sensititre ESBL antimicrobial susceptibility testing panel correctly differentiated all the ESBL-producing isolates from organisms containing other resistance mechanisms. The MIC was unaffected by inoculation directly from the selective *Brilliance* ESBL Agar compared to the non selective CBA.

Introduction

Extended spectrum β -lactamases (ESBL) are usually plasmid-mediated β -lactamases most commonly found in *Klebsiella pneumoniae*, *Escherichia coli* and other Gram-negative bacilli, and have emerged as a major threat worldwide with limited treatment options¹.

Rapid identification and confirmation of ESBL-producing organisms in the clinical setting is imperative in improving patient outcomes and preventing spread. This study was undertaken to determine if Thermo Scientific Sensititre[™] ESBL antimicrobial susceptibility testing panel (Thermo Fisher Scientific, figure 1) could be used with Oxoid *Brilliance*[™] ESBL Agar (Thermo Fisher Scientific, figure 2) for the rapid identification, confirmation and Minimum Inhibitory Concentration (MIC) -determination of suspect organisms

FIGURE 1. Sensititre antimicrobial susceptibility testing panels



Methods

Thirty one *K. pneumoniae* (21 ESBL-producers), 22 *E. coli* (14 ESBL-producers) and 13 other Enterobacteriaceae (either ESBL-producers or AmpC-producers), were inoculated onto *Brilliance* ESBL Agar and Columbia Blood Agar (CBA). All plates were incubated at 36±1°C for 24 hr.; any *Brilliance* ESBL Agar showing no growth at this time were re-incubated for a further 24 hr.

Sensititre ESBL antimicrobial susceptibility testing panel consists of the following antibiotic combinations: cefotaxime/clavulanic acid, ceftazidime/clavulanic acid, cefpodoxime/clavulanic acid, ceftriaxone/clavulanic acid and cefepime/clavulanic acid. Sensititre ESBL antimicrobial susceptibility testing panel inocula were prepared from colonies grown on both *Brilliance* ESBL Agar and CBA plates. QC testing was performed following the manufacturer's instructions and CLSI M100-S21².


Results

Brilliance ESBL Agar has previously proven to be a valuable screening method for detection of ESBL-producing Enterobacteriaceae, with a sensitivity and specificity of ≥95%³.

The Sensititre ESBL antimicrobial susceptibility testing panel correctly differentiated all ESBL-producing isolates that grew on *Brilliance* ESBL Agar from those containing other resistance mechanisms.

For 500 observations, the MIC was unaffected by inoculation of Sensititre ESBL antimicrobial susceptibility testing panel directly from selective *Brilliance* ESBL Agar compared to the non-selective CBA. The majority of observed MICs from these two procedures were within one doubling dilution, producing a 97% essential agreement.

FIGURE 2. *Brilliance* ESBL Agar



Conclusions

The Sensititre ESBL antimicrobial susceptibility testing panel can be used as a direct confirmatory test for suspect organisms isolated on *Brilliance* ESBL Agar. The prior use of *Brilliance* ESBL Agar does not influence the MIC result. The combination of *Brilliance* ESBL Agar with the Sensititre ESBL antimicrobial susceptibility testing panel offers a rapid, simple solution for identification, confirmation and antimicrobial susceptibility testing of ESBL-producing organisms.

References

1. Goyal, A., Prasad, K.N., Prasad, A., Gupta, S., Ghoshal, U., Ayvagari, A. (2009). Extended spectrum beta-lactamases in *Escherichia coli* and *Klebsiella pneumoniae* and associated risk factors. *Indian J Med Res.* 129:695-700.
2. Clinical Laboratory Standards Institute. Performance Standards for Antimicrobial Susceptibility Testing, Twenty-first informational supplement (M100-S21).
3. Huang, T-D., Bogaerts, P., Berhin, C., Guisset, A., Glupczynski, Y. Evaluation of *Brilliance*[™] ESBL Agar, a Novel Chromogenic Medium for the Detection of Extended-Spectrum Beta-Lactamases-Producing *Enterobacteriaceae*. *J. Clin. Micro.* 48:6 2091-2096;

© Oxoid Ltd. Photographs may not be extracted or reproduced in any way. Folio no. LT2036/March/2012.
This information is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others.