

---

**Cyto-Tek® 2500**  
Cytocentrifuge

*Imagine... consistently  
superior cytology  
preparation*



continuous innovation for pathology



## Cyto-Tek® 2500 Cytocentrifuge

The Cyto-Tek® 2500 Centrifuge provides optimal cell recovery and advanced preservation of cellular structure in non-gynecological monolayer slide preparation.

The 2500 was designed with enhanced features focusing on user safety, including the introduction of a new sealable rotor housing with lid which significantly reduces aerosol exposure.

The housing is completely removable and transportable, providing added convenience.

Technicians now have the ability to load and unload specimens in vented areas away from the centrifuge, under a fume hood.


The Cyto-Tek® 2500 Centrifuge has a sleek and modern design and provides the superior clinical performance users have come to expect from Sakura's Cyto-Tek® Centrifuges.

### Dimensions

W x D x H	375 x 490 x 225 mm
Weight	11.7 kg

### Features

Capacity	12
RPM	200 ~ 2.500 by 10rpm increments
Program memory	30
Autoclavable	Yes (Guard bowl)



## Items & options

4300	Cyto-Tek® 2500, 100-240V; 50/60 Hz	4333	Cyto-Tek® 1 mL Filter Paper; 200/pcs
4301	Chamber Kit 1mL; 50/pcs	4334	Cyto-Tek® 6/12 mL Chamber Cap; 25/pcs
4326	Cyto-Tek® 6/12 mL Base Holder; 12/pcs	4335	Cyto-Tek® 1 mL Chamber Cap; 200/pcs
4327	Cyto-Tek® 6 mL Gasket; 25/pcs	4336	Cyto-Tek® 1 mL Base Holder; 200/pcs
4328	Cyto-Tek® 12 mL Fluid Chamber; 12/pcs	4337	Cyto-Tek® 12 mL Gasket; 25/pcs
4329	Cyto-Tek® 1 mL Fluid Chamber; 200/pcs	4338	Cyto-Tek® 6 mL Filter Paper; 500/pcs
4331	Cyto-Tek® 6 mL Fluid Chamber; 12/pcs	4339	Cyto-Tek® 12 mL Filter Paper; 500/pcs

Visit our website [sakura.eu](http://sakura.eu)

Sakura Finetek Europe B.V., Flemingweg 10a, 2408 AV Alphen aan den Rijn  
P.O. Box 362, 2400 AJ Alphen aan den Rijn, The Netherlands

